



NEWLY LAUNCHED

Bio Vision

Vision of future and Optimum performance

Every Drop of Precision

- Easy PUSH & FLIP system
- Ergonomic design & Smooth Plunger
- Strong Grip with window on hand grip
- Clear digital display with lock
- Light weight - 70 gms.



BIOSYSTEM COMPANY

" Take Pipetting to Generation Next "



Bio Vision

Biosystem newly introduced a (Bio Vision) Model at very affordable price for pathology, medical college, professional educational, laboratories



Additional features of Classic Plus

- Comfortable advance finger rest
- Push button with soft rubber
- The Weight just 70 gm (Approximately)
- Instant volume setting
- Simple in lab calibration
- Economical price


Accessories



BIOTIPS

- Made from virgin plastic
- Fully autoclavable
- Universal microtips

Cat No.	Description	Pack / Size
BV/003/WS	White (0.2-10 μ l)	1x1000 pcs.
BV/003/Y	Yellow (2-200 μ l)	1x1000 pcs.
BV/003/B	Blue (200-1000 μ l)	1x1000 pcs.
BV/003/WM	White (0.5-5 ml)	1x100 pcs.
BV/003/WL	White (1-10 ml)	1x100 pcs.




BIOTIPS BOX

- Made from heavy acrylic
- Universal pipette stand fitting upto 4 pipette & with tips sprae.

Cat No.	Description	Pack / Size
BV/017/W	White (0.2-10 μ l)	1x96 pcs.
BV/017/Y	Yellow (2-200 μ l)	1x96 pcs.
BV/017/B	Blue (200-1000 μ l)	1x96 pcs.

BIOSTAND



Cat No.	Description	Pack / Size
BV/004/A	BIOSTAND 4holes	1 pcs.
BV/004/B	BIOSTAND 5holes	1 pcs.
BV/004/C	BIOSTAND 2holes	1 pcs.
BV/004/D	BIOSTAND 1holes	1 pcs.

Order Specification

FIXED VOLUME

Cat. No.	Volume	Inaccuracy	Imprecision	Tips to be used
BV/005/5	5 μ l	\pm 2.5	< 1.5	Yellow
BV/005/10	10 μ l	\pm 1.2	< 0.8	Yellow
BV/005/20	20 μ l	\pm 1.0	< 0.5	Yellow
BV/005/25	25 μ l	\pm 1.0	< 0.4	Yellow
BV/005/50	50 μ l	\pm 1.0	< 0.4	Yellow
BV/005/100	100 μ l	\pm 0.8	< 0.3	Yellow
BV/005/200	200 μ l	\pm 0.8	< 0.3	Yellow
BV/005/250	250 μ l	\pm 0.8	< 0.3	Blue
BV/005/500	500 μ l	\pm 0.8	< 0.3	Blue
BV/005/1000	1000 μ l	\pm 0.8	< 0.3	Blue
BV/005/5000	5000 μ l	\pm 0.8	< 0.3	White
BV/005/10000	10000 μ l	\pm 0.8	< 0.3	White

VARIABLE VOLUME

Cat. No.	Volume	Increment	Inaccuracy	Imprecision	Tips to be used
BV/006/200	0.5 - 10	0.1	\pm 5.0 - 2.0	< 5.0 - 1.0	White
BV/006/210	10 - 100	1.0	\pm 2.0 - 1.5	< 1.0 - 0.6	Yellow
BV/006/220	5 - 50	1.0	\pm 2.5 - 1.4	< 1.5 - 0.6	Yellow
BV/006/230	2 - 20	0.1	\pm 4.0 - 2.0	< 2.0 - 0.8	White
BV/006/240	20 - 200	1.0	\pm 2.0 - 1.5	< 0.8 - 0.4	Yellow
BV/006/250	100 - 1000	10.0	\pm 1.5 - 1.2	< 0.6 - 0.4	Blue
BV/006/260	0.5 - 5 ml	100.0	\pm 1.5 - 0.5	< 0.6 - 0.4	White
BV/006/270	1 - 10 ml	100.0	\pm 1.5 - 0.5	< 0.6 - 0.4	White

MULTI-CHANNEL VARIABLE VOLUME

Cat. No.	Volume	Inaccuracy	Imprecision	Tips to be used
BV/006/MV/50	5 - 50 μ l	\pm 5.00 - 1.50	< 2.00 - 0.80	Yellow
BV/006/MV/300	30 - 300 μ l	\pm 4.70 - 1.20	< 1.80 - 0.70	Yellow

NOTE : Special volume pipettes other than specified above may be supplied on demand.

BIOSYSTEM COMPANY

Office-Work: 105, Sahakar Industrial Estate, Navghar Road, Vasai (E), Thane 401 210, Maharashtra, India.

Tel: +91 250 2393979 / 3209901

Email: info@micropipette.in / biosystem2007@yahoo.com

Website: www.micropipette.in